




A new conception of Temporo Mandibular Joint (TMJ) and the method to assemble it

Technology



Benefits

- The glenoid insert can be easily fixed and release from the glenoid support, by a manual action on a tab, in order to facilitate it's precise adjustment or it's change to perfectly match with the condylar part.

Applications

- Prosthesis placed by maxillofacial surgeons in case of Temporomandibular joint dysfunction (TMJD)

Keywords

TMJ
 TMJD
 Temporo Mandibular Joint
 Metallic taste in mouth
 Maxillofacial surgery

Intellectual Property

WO 2018/149849
 published 08/23/2018
 EU, US, CA, CN, JP, BR

Development Status

Technical drawing and invention patented

Partnership

.Licencing or collaboration

contact

François-Xavier DENIMAL

Business Developer

+33 6 13 84 36 28

francois-xavier.denimal@sattnord.fr

find other technologies on
www.sattnord.fr



SATT Nord

25, avenue Charles St Venant – 59800 LILLE – France

+33 3 28 36 04 68 – tech@sattnord.fr