

Biaxial tests on granular materials and powders

BIAXIAL TESTS - GRANULAR MATERIALS - DEFORMATION

COMPETITIVE ADVANTAGES

- · Biaxial deformations;
- High precision on measurement displacements (of the order of 1/10e of μm);
- Possibility to perform tests on granular materials with fluids.

APPLICATIONS/MARKETS

- Study of granular materials and powders;
- Natural materials (sand, gravels, clay, corn, coffee): geotechnics, soil mechanics;
- Materials derived from industrial processes (pharmaceutical tablets, steel and glass beads, food pellets): pharmaceutical and food industry.

INTELLECTUAL PROPERTY

- Patent FR1161654;
- Collaboration sought: License for industrialization;
- Readiness level: laboratory prototype.

LABORATORY

 Laboratoire d'Etude des Microstructures et de Mécanique des Matériaux (LEM3)

CONTACT

Ludovic GOBY
Development officer
Materials, Processes, Chemistry
Tel.: 03.80.40.34.97 - 06.43.65.51.20
Mail: ludovic.goby@sattge.fr



PRESENTATION

The LEM3 laboratory has developed a device to test the behavior of granular materials and powders. It allows studying grain deformation and a granular assembly submitted to an external compression.

The device is made of a main box, which receives the tested material, on which are fixed longitudinal and lateral walls, partly mobile to apply pressure along those two axes. It also comprises a flexible membrane, on which pressure can be applied to increase the applied force intensity.

